

BURMEIER

Nursing bed "Arminia II"

Manufactured to:

- ▶ DIN EN Standard 1970:2000
- ▶ DIN EN Standard 60601-2-38/A1:2001
(limited to the requirements of domestic care)



CE

Model 51.3431.33

Nursing bed "Arminia II" – model 51.3431.33

► Design

For mattress size 90 x 200 cm, four-part mattress frame, mattress height adjusts from 40-80 cm, back and upper-leg support as well as mattress height adjustable by electric motors. Lowered-foot position can be altered by swinging down mattress frame. Metal parts sealed by polyester powder coating, colour: topaz, beech effect/solid beech wooden surround.

► Dimensions

Length 220 cm, width 102 cm

► Order numbers

51.3431.33 Nursing bed, colour: topaz

19.1400.1 Sit-up facility (included as standard)

Subject to technical modification without notice

► Accessory*



Serving tray

► Accessory*



Add-on safety frame rail

► Accessory*



Storage/transport aids

► Accessory*



Foam leather cover

► Accessory*



Reading lamps

► Accessory*



Bed extension, including safety frame rails (can also be retrofitted)

► Detail



Lower-leg rest lowered manually

► Detail



Nursing bed dismantled

► Detail



Lowered-foot position

* Accessories available at extra charge.