

TECHNICAL DATA

Overall length 207cm
 Overall width 100cm
 Sleeping surface 88 x 191 cm
 Spring height:
 Hi-Position 74 cm
 Low-Position 43.5 cm
 Hi-Low Travel 31 cm
 Back rest inclination 70°
 Thigh rest inclination 40°
 Castors diameters Ø100 mm
 Ventilated EG-sheet steel mattress platform
 Manual operated Hi-Low system
 Manual operated backrest
 Manual operated knee section
 Safety side rails
 Revolving wall bumpers



STANDARD EQUIPMENT
FOLDING SIDE RAILS


OPTIONS & ACCESSORIES

- I.V rod adjustable
- Patient chart holder
- Patient card holder
- Polyether foam mattress




RETRACTABLE CRANKING HANDLE


BED POSITIONS




Backrest



Fowler



Low-Position



Hi-Position

Rehab Mart



Quality Healthcare Furniture



Quality Healthcare Furniture



REHABILITATION BED

B-1.3