



Posey Abduction Wedges

POSEY ABDUCTION WEDGES ARE LIGHTWEIGHT AND PROVIDE FULL LEG AND ABOVE THE KNEE ABDUCTION. THE UNIQUE CONSTRUCTION OF AN OUTER VINYL MESH AND A BREATHABLE FABRIC WITH A POLYSTYRENE BEAD FILL ALLOWS AIR TO FLOW THROUGH THE WEDGE, ELIMINATING THE DISCOMFORT OF CONVENTIONAL FOAM FILLED PRODUCTS, WHICH BUILD UP HEAT. BOTH WEDGES FEATURE A ZIPPERED CLOSURE, ENABLING BEADS TO BE ADDED OR REMOVED FOR PATIENT ACCOMMODATION. 4" WIDE SOFT FABRIC STRAPS WITH HOOK CLOSURES KEEP THE WEDGES COMFORTABLY SECURED BETWEEN THE KNEES. THE SMALL ABOVE KNEE WEDGE (CAT. NO. 6302S) ALLOWS COMFORTABLE ABDUCTION OF THE THIGHS WHILE SEATED. THE LARGE FULL LEG WEDGE (CAT. NO. 6302L) SECURES THE LEGS AND COMES WITH TWO SOFT FABRIC STRAPS WITH HOOK FASTENERS FOR EXTRA STABILITY. MACHINE WASHABLE. 6302L 22"L X 5"H AND TAPERS FROM 15" AT THE FEET TO 7" AT THE THIGHS