



Posey Elevation Wedges

POSEY ELEVATION WEDGES OFFER COMFORT AND SUPPORT FOR PATIENT'S REQUIRING ELEVATION OF THEIR TORSO, LEG AND/OR ARM. THE LIGHTWEIGHT WEDGES ARE CONSTRUCTED OF AN OUTER VINYL MESH AND A BREATHABLE FABRIC FILLED WITH POLYSTYRENE BEAD FILLS. THE UNIQUE DESIGN ALLOWS AIR TO FLOW THROUGH THE WEDGE, ELIMINATING THE DISCOMFORT OF CONVENTIONAL FOAM FILLED PRODUCTS, WHICH BUILD UP HEAT. ALL THREE WEDGES FEATURE A ZIPPERED CLOSURE, ENABLING BEADS TO BE ADDED OR REMOVED FOR PATIENT ACCOMMODATION. THE ELEVATION TORSO WEDGES (CAT. NO. 6311L) ARE IDEAL FOR USE WITH ALTERNATING PRESSURE MATTRESSES AND HELP PREVENT CLIENT SLIPPAGE. THE ELEVATION ARM WEDGE (CAT. NO. 6311A) ELEVATES THE FOREARM FOR HAND AND WRIST SUPPORT AND IS IDEAL FOR USE WITH CVA PATIENTS, OR THOSE ON AN IV. MACHINE WASHABLE. 6311 ELEVATION TORSO WEDGE- REG. (35/55 DEGREE ANGLE) 18"L X 18"W X 13"H 6311A ELEVATION ARM WEDGE- (30/60 DEGREE ANGLE) 18"L X 11"W X 7 1/2 "H