

Posey Lateral Wedges



POSEY LATERAL WEDGES ARE LIGHTWEIGHT AND PROVIDE BOTH POSITIONING SUPPORT AND COMFORT. WEDGES CAN BE EASILY PLACED IN POSITION UNLIKE FOAM PILLOWS, AND IT IS POSSIBLE FOR A SINGLE CAREGIVER TO TURN THE PATIENT WITHOUT ASSISTANCE BY USING THE WEDGES. THE LATERAL WEDGES ARE IDEAL FOR KEEPING-BED RIDDEN PATIENTS REPOSITIONED, ELIMINATING PRESSURE POINTS OVER THE BONY PROMINENCES AND REDUCING THE INCIDENCE OF PRESSURE ULCERS. THE UNIQUE CONSTRUCTION OF AN OUTER VINYL MESH, BREATHABLE FABRIC SHELL AND A NON-ABSORBENT, POLYSTYRENE BEAD FILL ALLOWS AIR TO FLOW THROUGH THE WEDGE, ELIMINATING THE DISCOMFORT OF CONVENTIONAL FOAM PILLOWS, WHICH BUILD UP HEAT. ALL THE WEDGES FEATURE A ZIPPERED CLOSURE, ENABLING BEADS TO BE ADDED OR REMOVED FOR PATIENT ACCOMMODATION. MACHINE WASHABLE. 6309M PATIENTS UNDER 180LBS. 16"L X 6"W X 7"H