



HE7720 ANTI-DECUBITIS AIR PUMP

The Wind Star series could separately go with following different pads: bubble mattress (regular/micro air holes), as a wide of mattress selection to provide patient can select comfortable mattress to release heal bedsore pressure.

Pump specification:

HE7720002V Adjustable Pressure Control

1. Size: 27 x 14 x 13 cm
2. Weight: 1.6KG
3. Power: 115V/220V
4. Pressure range: 110-140mmHg
5. Air output: 3.2L
6. Cycle time: every 5 minutes (bubble mattress)
7. Power cord length:1.8M

HE7310 (MICRO AIR HOLES) BUBBLE MATTRESS

The mattress are designed to match with pump's air output-input cycle and formation mttress air channel up and down in every 5 minutes to makes air can go through patient's back to release pressure sore occur.

As this idea, the mattress effective match with air pump operate can reduces the workload for the caregiver and provide patient more comfortable, helping to prevent heal bedsores.

Based on several situation to use, the mattress we are separate in 3 kinds of different model mattress for user's to select and each model has 2 styles (regular/micro air holes) for user's has more wider range to use.

The micro air holes can via air loss to provide patient's more dry environment to prevent patient's Back due to long term lying on the bed and sweat can't evaporation caused the bed soared occur BUBBLE MATTRESS SYSTEM.

It is an economical type and can be directly used on the hospital bed or nursing home, also easy user to put bubble mattress on the other mattress.

This mattress are suitable with our air pump: HE7500/HE7510 series, HE7600/HE7610 series, and HE7700/HE7710 series.

1. MATERIAL: PVC(EN-71, NON-TOXIC,
2. COLOD RESISTANT:-25 DEGREE)
3. THICKNESS:0.35mm
4. SIZE:78" x 34"(INFLATTED)
5. HEIGHT OF CELLS (INFLATED): 2 1/2"
6. NUMBER OF CELLS (INFLATED): 130
7. MICRO AIR HOLES: 46 CELLS

