



PH649 Permaflow Fibre Filled

PL020 Permalux Fibre Filled

Overlay Mattress

Description :

- A fibre filled Overlay mattress
- PL020-with integral stitched Permalux waterproof cover

User Weights :

- Minimum - 6 Kgs
- Maximum - 95 Kgs

Dimensions:

- Length - 195 cm
- Width - 90 cm
- Depth - 15 cm
- Weight - 7.5 Kgs

Mattress Cover (PL020):

Material

- Polyurethane coated fabric
- Impregnated with *Acti-fresh* biocide
- 2 way stretch fabric
- Vapour permeable
- Waterproof
- Flame resistant – (Crib test 5)

Cleaning

- Hand clean using a 1000 p.p.m. solution of Sodium Hypochlorite - Rinse well
- Machine washable
- Maximum temperature 71 degrees Celcius

Warranty :

- 1 Year manufacturers warranty.

Due to its policy of continued development and improvement, Park House Healthcare Ltd reserves the right to amend any details and specifications without notice.