

Karma Tilt-in-Space Children Wheelchair CP-33



The design and development of a personal CP wheelchair is quite a challenge for every wheelchair manufacturer. The all-new KM-CP30 is the proud result of three years of collective efforts between Karma and senior German medical engineers. Continuing the great traditions of its predecessors, the 3rd generation Karma CP chair boasts many revolutionary design features that are second to none.

- Fixed seat and backrest system; ready for later upgrade
- Tilt-in-space seat offers pressure relief and a shift in center of gravity
- Modular design offers quick and convenient assembly, folding and disassembly
- Compact when folded, requiring minimal storage space
- Frame weights only 11kg once folded; nimble and easy to transport



Dining Tray

- Can be fine-tuned from 0° to 45° angles for enhanced reading and writing comfort



Separator

- Adjustable separator can separate thighs from further damage



trunk support

- Adjustable trunk support protects body from leaning or distortion




Safety belt/

- Multi-functional safety belt provides added security and protection for users

Specifications

Unit : c m

*Size: 9"-16 " *F/R Wheels: 5"-17"/14"-22" *Max. User Weight: 80kg



23-41



23-46




40-57



43.5




15-30.5



83-106



48-48.5



by customer



by customer



11kg/20kg
± 1cm ± 0.5kg