

Buggy



Options and Accessories



Padded head support



Sunshade



Trunk support



Pad on side rail



Thigh support

HE3800061U

- * For children with special needs
- * Lightweight and easy to fold
- * Features adjustable seat depth, seat width and angle adjustable footplates. Allows the buggy to grow with the child
- * Provides maximum comfort and flexibility

Specification

- * Rear seat-to-floor heights (Adjustable): 8" wheel, 12"-14"
- * Front seat-to-floor heights (Caster size): 6", 18"-20"
- * Knee-to-heel range (Front seat-to-footplate adjustable): 8"-14"
- * Depth of seat (Adjustable): 12"-16" (Range adjusts in 2" increments)
- * Width of seat 14.5"
- * Footplate angle (Adjustable): 0°-45° (Range adjusts in 12° increments)

Rehab Mart

www.RehabMart.sg T:62500555