



Maclaren Major Elite Buggy

Maclaren Major Elite special needs pushchair is considered THE everyday buggy for children with special needs. In most respects, it looks and performs like a regular buggy with swivel front wheels and square-grip handles, but it incorporates an extra-strength, square-section aluminium chassis making it super-sturdy yet lightweight. The brakes on each rear wheel have positive-lock safety catches, and a five point safety harness come as standard.

The three-position detachable, rigid foot rest ensures the Major Elite can grow with your child into adulthood if required - as it has been tested by the Medical Devices Agency up to 50 Kg. With this kind of carrying capacity as well as a hard-wearing, washable seat cover, the Major Elite will withstand years of rugged service. Dimensions : Seat Width 15"(38cm), Seat Depth 10"(25cm), Seat Height 20"(51cm) , Adjustable Leg Height 9"/11"/15.5"(23/28/39cm) It weighs 6.75kg, (7.5kg with the footplate) and folds easily to to 46" x 9 3/4" x 8 1/4" (117cm x 25cm x 21cm).

Optional Accessories -

Shopping bag, Matching sun shade, Transparent rain cover (must be used with sunshade), Chest pad (to prevent older children undoing the safety harness), Two padded lateral inserts (for extra support which can be fitted left and right of the harness) A padded seat cover for additional comfort, A padded footfuff for warmth and weather protection

Rehab Mart