

# CHILDREN'S WHEELCHAIR X2 JUNIOR (MODEL 2.351-351)

HMV-number: 18.50.03.1050



**MEYRA**  
ORTOPEDIA

Adjustment	Operation	Accessories
<ul style="list-style-type: none"><li>• Extremely light and compact folding wheelchair with slanted front chassis</li><li>• Active provision concept with adjustment quality from seat width 28–46 cm (11.0-18.1 inch) and seat depth 32-38 cm (12.6-15.0 inch), lower shank length from 30 cm (11.8 inch)</li><li>• Many adjustment possibilities to functional changes</li><li>• especially light running through fixed frame driving behaviour</li><li>• Wheel camber adjustable</li></ul>	<ul style="list-style-type: none"><li>• Variable back angle adjustment up to 15°</li><li>• Push handles optionally height adjustable</li><li>• Arm supports removable, lockable and swivel-away, optionally height adjustable</li><li>• Continuous footboard or single footplates depth adjustable</li><li>• Adjustable lower shank length</li></ul>	<ul style="list-style-type: none"><li>• Brake lever extension</li><li>• Heel cup</li><li>• Special air pump</li><li>• X-trends accessories</li><li>• Accessory net</li></ul>

## FIELDS OF APPLICATION

### Child wheelchair X2 Junior model 2.351-351

#### Fields Of Application

- For active larger kids or adolescents with the provision goals enhancement, compensation, independence and mobility through exact, demand oriented adaptation
- For long term and permanent provision at home, in school and during recreation

## DETAILS



#### Everything for the feet

With single footplates, automatically folding footboard or the continuous footboard, for adolescents or adults, with footboard locking device or without - the X 2 junior is master of any demand!

## X2 JUNIOR DIMENSIONS

### Technical Data

Seat width	28 / 11	30 / 11.8	32 / 12.6	34 / 13.4	36 / 14.2
Width unfolded <sup>1)</sup>	45 / 17.7	47 / 18.5	49 / 19.3	51 / 20.1	53 / 20.9
Width folded	30 / 11.8 (0.5° wheel camber)	30 / 11.8 (0.5° wheel camber)	30 / 11.8 (0.5° wheel camber)	30 / 11.8 (0.5° wheel camber)	30 / 11.8 (0.5° wheel camber)
Seat depth	32 / 34 / 36 / 38 // 12.6 / 13	32 / 34 / 36 / 38 // 12.6 / 13	32 / 34 / 36 / 38 // 12.6 / 13	32 / 34 / 36 / 38 // 12.6 / 13	32 / 34 / 36 / 38 // 12.6 / 13
Seat height front/rear	44 to 55.5 / 40.5 to 52 // 17.	44 to 55.5 / 40.5 to 52 // 17.	44 to 55.5 / 40.5 to 52 // 17.	44 to 55.5 / 40.5 to 52 // 17.	44 to 55.5 / 40.5 to 52 // 17.
Overall height	80 / 85 / 90 / 95 // 31.5 / 33	80 / 85 / 90 / 95 // 31.5 / 33	80 / 85 / 90 / 95 // 31.5 / 33	80 / 85 / 90 / 95 // 31.5 / 33	80 / 85 / 90 / 95 // 31.5 / 33
Back height	variable from 30 to 47.5 / 11.	variable from 30 to 47.5 / 11.	variable from 30 to 47.5 / 11.	variable from 30 to 47.5 / 11.	variable from 30 to 47.5 / 11.
Arm support height	19 - 25 / 7.5 - 9.9	19 - 25 / 7.5 - 9.9	19 - 25 / 7.5 - 9.9	19 - 25 / 7.5 - 9.9	19 - 25 / 7.5 - 9.9
Length with foot support	81 / 31.9	81 / 31.9	81 / 31.9	81 / 31.9	81 / 31.9
User weight in kg / pd	120 / 264.5	120 / 264.5	120 / 264.5	-120 / 264.5	120 / 264.5
Empty weight	11 / 24.2	11 / 24.2	11 / 24.2	11 / 24.2	11 / 24.2
Permissible overall weight	131 / 288.7	131 / 288.7	131 / 288.7	131 / 288.7	131 / 288.7
Transport weight	7.5 / 16.5	7.5 / 16.5	7.5 / 16.5	7.5 / 16.5	7.5 / 16.5
Driving wheel	22" / 24"	22" / 24"	22" / 24"	22" / 24"	22" / 24"

Seat width	38 / 15.0	40 / 15.8	43 / 16.9	46 / 18.1
Width unfolded <sup>1)</sup>	55 / 21.7	57 / 22.5	60 / 23.6	63 / 24.8
Width folded	30 / 11.8 (0.5° wheel camber)	30 / 11.8 (0.5° wheel camber)	30 / 11.8 (0.5° wheel camber)	30 / 11.8 (0.5° wheel camber)
Seat depth	32 / 34 / 36 / 38 // 12.6 / 13	32 / 34 / 36 / 38 // 12.6 / 13	32 / 34 / 36 / 38 // 12.6 / 13	32 / 34 / 36 / 38 // 12.6 / 13
Seat height front/rear	44 to 55.5 / 40.5 to 52 // 17.	44 to 55.5 / 40.5 to 52 // 17.	44 to 55.5 / 40.5 to 52 // 17.	44 to 55.5 / 40.5 to 52 // 17.
Overall height	80 / 85 / 90 / 95 // 31.5 / 33	80 / 85 / 90 / 95 // 31.5 / 33	80 / 85 / 90 / 95 // 31.5 / 33	80 / 85 / 90 / 95 // 31.5 / 33
Back height	variable from 30 to 47.5 / 11.	variable from 30 to 47.5 / 11.	variable from 30 to 47.5 / 11.	variable from 30 to 47.5 / 11.
Arm support height	19 - 25 / 7.5 - 9.9	19 - 25 / 7.5 - 9.9	19 - 25 / 7.5 - 9.9	19 - 25 / 7.5 - 9.9
Length with foot support	81 / 31.9	81 / 31.9	81 / 31.9	81 / 31.9
User weight in kg / pd	120 / 264.5	120 / 264.5	120 / 264.5	120 / 264.5
Empty weight	11 / 24.2	11 / 24.2	11 / 24.2	11 / 24.2
Permissible overall weight	131 / 288.7	131 / 288.7	131 / 288.7	131 / 288.7
Transport weight	7.5 / 16.5	7.5 / 16.5	7.5 / 16.5	7.5 / 16.5
Driving wheel	22" / 24"	22" / 24"	22" / 24"	22" / 24"

• Dimensions in cm / inch (1:2,54) • Weights in kg / pd (1:0,45) • Tolerance ± 1,5 cm / 0,59 inch • Speed in kmh / mph (1:1,61) •