



# Rehab Mart



## HAND 701

KM-7501 exclusively for children

### Caring for children's world

Carer-friendly. All-around protection.  
Maneuver flexibility.

#### HAND 701,

Handle your children with great care.



#### Here-we-go handle

Antelope horn-shaped handle is ergonomically designed for better gripping, which makes it easy to steer and push the chair. Its lively look also represents the quality of children.



#### Outward extended front wheels

Front wheels are extended outwardly to significantly increase contact surface with the ground, making its body more stable and much safer in propelling.



#### Caring footrest

Inward curving footrest is intelligently designed, not only stabilizes lower body positioning, but provides greater protection by avoiding direct collision and enhances rotation flexibility as well.



#### Vertical footrest

The vertical footrest is ergonomically designed to adjoin children's body type, allow their legs to be placed in the correct position, and make sitting more relaxed and cozy.



#### One-piece footplate

The one-piece footplate with animated figures has greatly improved stability; combined with its easy-folding design, it provides children with joy and safety at the same time.



#### Ultra lightweight and compact

Weighting only 9.3 kg, the chair folds easily and is extremely compact after folding, which makes it easy for storing and traveling.

### Recommendation



#### Side-lifting transparent tray

The tray itself can be stored next to the armrest, thus greatly improves children's maneuvering ability. And its transparent design allows the carers to see the ground clearly, prevent accidents from happening, and enhance children's safety.

### Optional accessories

- Anti-tippers, provide greater sitting safety.
- Quick release 20" rear wheels, assembly with great convenience.
- 20" rear wheel cover, luminous designed and prevents hand pinching, greatly improves safety, even for self-propelling.



Width	F/R wheels												Max. User Weight
11"	6"-14F"	28	30	45	34	18	70	39	36	102	9.3	60kg	
13.5"	6"-14F"	34	30	51	34	18	70	39	36	102	9.3	60kg	

±1cm ±0.5kg

