

recline & elevating footrest come standard!



KM-5000F14

**KM-5000**

Width	F/R Wheels												Max. User Weight
17"	8"-24"	43cm	42cm	64cm	44cm	22cm	122cm	51cm	45cm	93/129cm	20kg	115kg (255lbs)	± 1cm ± 0.5kg

KARMA reserves the right to modify information herein without further notice.

**[Features]**

- Ultra-lightweight aircraft-grade aluminum alloy frame
- Standard seat width: 17"
- Weight capacity: 115kg
- Detachable armrests and elevating footrests
- Foldable frame with anti-tippers
- Backrest reclines smoothly from 90°-163°
- Detachable and height-adjustable headrest
- Product weight: 20kg

**[Options]**

- Hub brake
- 24" flat free rear wheels
- Detachable and washable cushion
- Armrest:
  - ① Flip-back, detachable, height adjustable armrest
  - ② Flip-back, detachable, armrest
- Footrest:
  - ① Amputee footrest
  - ② Swing-away footrest

