

KP-65 SERIES

The all-new KP65 series is Karma's 5th generation of power chairs. True to its predecessors, reputed for high-torque and reliability, the new KP65 series allows for easy upgrades and includes advanced features to compliment the users individual needs.



Special Features P65TR

- Adaptive frame for seating upgrades.
- Lightweight aircraft-grade aluminum alloy 6061.
- Multi-functional with power tilt and continuous adjustable power recline.
- Accommodates:
 - Deluxe Captain Seat.
 - Custom seating for varying posture and pressure needs.
 - Most after-market seating, such as Jay cushions and backs.
- High efficiency motor package with excellent torque and reduced running noise levels.
- Great maneuverability with added safety.
- Comfortable and pressure-relieving seating.
- Advanced VSI controller, adaptable for future power tilt and recline upgrades.

KP-65

- Deluxe Captain Seat
- Fixed seat
- Manual recline

KP-65T

- Deluxe Captain Seat
- Power tilt-in-space
- Manual recline



KP-65TR adapting with KARMA standard cushion and back

- available for adapting after market seating systems.
- Power tilt-in-space
- Power recline



For example,

KP-65TR adapting with jay cushion and back

- Power tilt-in-space
- Power recline



Backrest Design for Comfort (Recline)

- Accommodates various seating systems and offers either manual or power recline.
- Continuous recline offers comfortable positions ranging from 90° to 165°.
- Adjusts smoothly and safely to desired position, relieving pressure on seat and trunk.

Seat Design for Comfort (Tilt-in-Space)

- Accommodates various seating systems and offers power tilting.
- Easily and smoothly adjusts to the most desirable position, shifting user's center of gravity to maximize comfort.
- Adjustable seat depth to better fit user's leg length.

Free-wheel Hubs



Free-wheel hubs allow safe, mechanical handling; even if battery power is completely depleted, the chair can be manually pushed when the hub is gently disengaged.

Shock Absorbers

Independent rear wheel suspension and dual shock absorbers ensure a smooth ride over uneven terrain. A well balanced chassis design provides maximum comfort and minimum pressure even during extended use.



Separate base unit and seating systems provide therapists the freedom to prescribe the most beneficial seating systems.



Advanced VSI control system delivers safe transitions between driving and seating modes.



Multi adjustable footrest allows the user to choose the most comfortable leg position. Also available in power elevating options.



CE AS FDA



International quality certificates



Founded in 1987 Karma is the leading manufacturer aluminum-alloy wheelchairs in Taiwan and one of the leading manufacturers of lightweight power wheelchairs and scooters. An exporter of its products to 35 countries Karma works closely with its network of dealers to offer a simple and accurate policy. Simply put we strive to offer the best quality, most innovative product to our agents and their clients at an affordable price.

We have recently :



- Strategically relocated headquarters to a newly established state-of-the-art office building and manufacturing plant, thereby integrating R&D and manufacturing together.
- Expanded our international sales and marketing to include distribution of " **KARMA** ", " **SOMA** " and " **Discovery** " brand name products in 30 countries.
- Accumulated 30 patents globally.



Specifications

Description	Model	KP-65	KP-65T	KP-65TR
Standard seat width		18" or 20" (Deluxe Captain Seat System)		
Optional with extra charge		16"/18"/20"/22" (Interface Seat System)		
Frame		Aircraft grade aluminum alloy		
Batteries		12V 50 A/H X 2pcs		
Speed		1~13 km/hr		
Motor		DC 400W X 2pcs		
Uphill		15 degrees (for safety, do not exceed 12 degrees)		
Max range		45km (distance and operating time can vary according to conditions)		
Brake		Electromagnetic brake		
Controller		50 AMP		
Charger		Input : 115 - 230 V AC Output : 24 V DC, 5A / hr		
Wheel size		Front caster: 2.80/2.50-4(diameter:22cm) Rear wheel: 3.00-10(diameter:40cm)		
Tire pressure		Front caster: 50 psi Rear wheel: 36 psi		
Shock absorber		Dual shock absorber with five tension levels		
Shock system		Independent suspension system		
Power tilt-in-space		X	O	O
Recline	Manual	O	O	X
Functional	Power	X	X	O
Overall length		125cm (with swing-away footrests and anti-tip)	102cm (w/o)	
Overall width		64.5cm		
Overall height		125cm		
Turning radius		135cm (with swing-away footrests and anti-tip)	95cm (w/o)	
Weight	with batteries	93kg	97kg	101.5kg
	w/o batteries	64kg	68kg	72.5kg

±1cm ±0.5kg



Headquarter in Taiwan



Thailand/Soma manufacturing operation



Karma United Kingdom

NATIONAL AGENT :

Rehab Mart

www.RehabMart.com.sg

Balestier: 12 Boon Teck Rd, #01-02 Shaw Corner

T: 62500555

Outram Park MRT: 100 Eu Tong Sen St #02-34

Pearls Ctr T: 63247077

Paya Lebar MRT: 10 Eunos Rd 8 #B1-02 Sing

Post Ctr T:6455 1551

SGH: Blk 7 L1 T:6221 1110

Fax: 65-62507576 Email: rehab@pacific.net.sg

KARMA MEDICAL PRODUCTS CO., LTD.

Tel: 886-5-2066688 Fax: 886-5-2067788

Http://www.karma.com.tw E-mail: eukarma@ms29.hinet.net