

# NETTI III COMFORT

The comfort chair, which has tilt, recline and the option of angle adjustable foot supports.



The NETTI III has been designed for users needing a high level of comfort and pressure relief. Netti III is a good alternative if individual adjustment and adaptation is essential for good sitting. A unique combination of seat and back cushions, ergonomic design of the frame and numerous adjustment possibilities make the Netti III a complete comfort wheel chair.

The tilt and recline functions allow for positioning according to the needs and wishes of the user. The forward tilt is especially useful for more active positions such as approaching a table, or entering the chair.

The chair allows the possibility to rest in a reclined position. Simple and durable functions also make the chair easy to use for those caring for the user. The chair can easily be folded for transport in car.

The chair can (in seat width 40 – 50 cm) be fitted with electrical control for tilt, recline and foot support functions.

## ALU REHAB AS

Fabrikkeveien 22 Kvål  
NO-4323 Sandnes, Norway  
Phone: +47 51 68 18 88  
Fax: +47 51 68 18 89

## NETTI III COMFORT COMES WITH THE FOLLOWING EQUIPMENT AS STANDARD:

<p>→ <b>BACK</b></p> <p>→ <b>ARMRESTS</b></p> <p>→ <b>HEAD SUPPORT</b></p> <p>→ <b>BRAKES</b></p> <p>→ <b>FOOTRESTS</b></p> <p>→ <b>MAIN WHEELS</b></p> <p>→ <b>PUSH RIMS</b></p> <p>→ <b>FRONT CASTORS</b></p> <p>→ <b>ANTI TIP</b></p>	<p>Standard back with Velcro. The back rest can be regulated by the user or companion. The backrest can be folded down for transport.</p> <p>With skirt guards, height adjustable and removable. The armrest pads can be adjusted in depth.</p> <p>Support for the head is always necessary on a chair with tilt and recline. There are 3 different headsupports to choose from. All head supports can be adjusted in height, depth and angle.</p> <p>The brake is mounted on the top frame tube (adjustable).</p> <p>Angle adjustable, height adjustable, swingable and removable.</p> <p>24" with quick release (In seat width 55 cm and 60 cm it comes with 24").</p> <p>Black plastic.</p> <p>7" with quick release (In seat width 55 cm and 60 cm it comes with 7").</p> <p>Height adjustable and length adjustable. Can be turned around and used as step climber.</p>
--	--



The Netti III can be tilted backwards and forwards.



The configuration of the wheelchairs differs from market to market.  
For other accessory: See pages 34-42.

### EXPLANATION TO THE TABLE BELOW:

\* The seat depth is measured without back cushion from front of the seat plate to the back tube placed in standard position. Under the chair there is an additional locking head for the gas spring. On seat depth 35-40 cm a short one is used, and on seat depth 43-60 a long one is used.

\*\* The seat height is measured from floor to top front seat plate.

## TECHNICAL SPECIFICATIONS:

CHASSIS SIZE	FUNCTIONAL SEAT WIDTH	SEAT DEPTH* Standard (Adjustable in depth)	SEAT HEIGHT**		BACK HEIGHT (From seat plate to top Velcro)	WEIGHT KG
			Standard (Without change of parts)	(With change of parts; 20"-24")		
35	37,5	45,5 (43-53)	44,5(47)	39,5-47	48	29,9
38	40,5	45,5 (43-53)	44,5(47)	39,5-47	48	30,4
40	42,5	45,5 (43-53)	44,5(47)	39,5-47	48	30,8
43	45,5	45,5 (43-53)	44,5(47)	39,5-47	48	31,9
45	47,5	50,5 (43-53)	44,5(47)	39,5-47	48	32,2
50	52,5	50,5 (43-53)	44,5(47)	39,5-47	48	33,1
55	57,5	50,5 (43-53)	44,5(47)	39,5-47	48	37,7
60	62,5	50,5 (43-53)	44,5(47)	39,5-47	48	38,8



The Netti III can be delivered with a gliding seat which makes positioning of the user easier.



The Netti III can be delivered with electrical functions.



The Netti III with electrical solutions comes with an integrated battery charger.

