

# MMED6000RD - MD300C FINGER TIP PULSE OXIMETER

## Technical Specification

### Patient Range

Adult, Pediatric

### SpO<sub>2</sub>

Range 0%-99%

Resolution 1%

Accuracy 70%-100%: ±2%

0%-69%: unspecified

Method Dual wavelength LED

### Heart (Pulse) Rate

Range 30-240 bpm

Resolution 1bpm

Accuracy ±2bpm or ±2%

### Display

Type LED display

Parameters SpO<sub>2</sub>, PR, Pleth bar.

Information Low Battery

### Environmental

Operating Temperature 5°C-40°C

Storage Temperature -10°C - 40°C

Operating Humidity 15%RH - 95%RH



Storage Humidity 10%RH - 95%RH

### Classification

Type of protection internally powered equipment

Safety IEC Standard 61601-1-1

### Michanical

Size 58 x 32 x 34 mm

Weight 50g (including 2 batteries)

### Battery

Type 2 AAA size 1.5V Alkali Batteries

or

2 AAA 1.5V rechargeable batteries

Operation time: minimum 30 hours for normal Operation



DISTRIBUTOR:

*Rehab Mart*



**Beijing Choice**

**Electronic Technology Co., Ltd.**

Address: Bailangyuan Building B1126-1128, Fuxing Road A36,  
Beijing 100039 P.R.China

Tel: +86 10 88203520, 88203551, 88204630, 88204631,

Fax: +86 10 88204632, 88204188

E-mail: [michael\\_zhang@choicemmed.com](mailto:michael_zhang@choicemmed.com)

[www.choicemmed.com](http://www.choicemmed.com)