



# SuperVac medi-pump



The Merridian SuperVac medi-pump offers medium to high suction aspirations. Designed for a wide variety of medical suction applications, this portable aspirator is ideal for hospitals, nursing homes, homecare and clinics.

## **SPECIFICATIONS:**

- Type: Diaphragm Type
- Motor: 1/8 HP AC 220—240V 50 Hz
- Airflow: 38—42 lpm
- Vacuum: 580 mmHg
- Weight 3.9Kg
- Size: L33.5 X W16 X H21.5 cm
- Equipped with adjustable vacuum gauge