



L-3



F- Floor to Push-handle
 B-Seat to back height
 C-Arm height from Seat
 D-Seat to Leg rest
 A- Floor to Seat height



E(b)-Distance bet: Outer Arms
 E(a)-Distance bet: Inner Arms
 K(a)- Width of cut-out
 K(b)-Length of cut-out

Model: L-3

A: 46cm-52cm

D: -

F: 123cm-129cm

I: 42cm

B: 79cm

E (a): 54cm

G: 41cm

J: 50cm

C: 23cm

E (b): 66cm

H: -

K (a): 32cm

(b): 39cm