



Dyna Self-propel Light Weight Aluminium Alloy Commode DL10598

Frame : powder-coated aluminium alloy

Armrest : Fixed height, flip-up

Seat : Plastic, 19" Width. Aperture size : 9 ¾" X 7"

Backrest : Plastic. Height -15"

Base clearance : 17"

Rear Wheels : 24" with air tyres Front Castor : 6" solid, swivel.

Legrests : plastic, retractable

Push-Handle height : 37"

Overall Dimension : 26'Wx31'Dx16.5'H

Rehab Mart