



**HE10665092 Lifting Commode 4 Casters**

- \* Push the user up from a seated position
- \* Height adjustable from ground to backrest: 92cm
- \* The distance between right-and-left leg tubing: 51cm
- \* The distance between front-and-rear leg tubing: 67cm
- \* Material: PU cushion, PU backrest, PU armrest
- \* Pail, spreadguide included
- \* Weight: 11Kg



HE10665092

**HE10666000 4 in 1 Reclining Shampoo Chair**

- \* Ideal for shower, shampoo, commode and transportation
- Reclining backrest
- \* Height adjustable from ground to backrest: 106cm
- \* The distance between right-and-left leg tubing: 56cm
- \* The distance between front-and-rear leg tubing: 84cm
- \* Material: PU cushion, PU backrest, PU armrest
- \* Pail, spreadguide included
- \* Weight: 12.9Kg



HE10666000

Adjust Reclining Backrest

*Rehab Mart*

www.RehabMart.sg T:62500555