






Order No.		Customer		Purchase Order	
Model Number	KM-AT_____	Issued Date		Sales Rep.	
Order Quantity	_____Set(s)				
Color Selection	<b>Frame Color :</b> Rose    Bright Yellow    Metallic Silver    Pearl Blue <b>Upholstery Color :</b> Charcoal    Karma Black				
Specification	Standard Accessories	Specification/Size			
	Standard				
	Seat Width & Depth	13 " x14 "    13 " x15 "    14 " x14 " 14 " x15 "    15 " x16 "    16 " x16 " customized _____ (Additional charge)			
	Seat Height	Standard(19.5 " ) -1 " ,    +1 " Customized: _____ (Additional charge)			
	Backrest Height	Standard ( 10 " ) ,    +5cm,    +10cm, Customized_____ (Additional charge)			
	Rear Wheel Angle	Standard 0° Option    (Additional charge)    3°    5° Option    (Additional charge)    7°    9°			
	Calf Length	Standard Height 15 " -1 " ,    +1 " ,    +2 " Customized____ " (Additional charge)			
Please refer to KMAT Catalogue	Upholstery	Fixed Adjustable ( Additional charge )			
	Rear Tire	PVC 5 " Gray Tire PU 3 " or 4 " or 5 " Elastic tire (Additional charge) 4 " Flashing Tire (Additional charge) 6 " Pneumatic Tire (Additional charge)			
	Footplate	Swivel Footplate Joined Footplate			
	Rear Wheel Set	Spoke Wheel Magnesium Alloy Wheel ( Additional charge ) Color:    Silver /    Black			

	Armrest Tube		W/O Armrest With Armrest (Additional charge)
	Rear Tire 24 " *1 "		Grey Color Color Tire (Additional charge) Color : Bright Yellow

Remarks : Standard parts - air pump, tool bag, anti-tips, push-handles. \_\_\_\_\_

Customer's Signature \_\_\_\_\_ Date \_\_\_\_\_