

Order Form

Purchase order Quotation

Address:
Shipping address:
Customer/ Remarks:

Active Wheelchair

3.350 X1



admissible user weight 120 kg

Price of basic model 3.350 (without armrests, legrests, etc.)

Sing Dollar conversion : **1.928,00 €**

Exampel: 3.350 - 40 - 54 - 101 - 183 - 287 = 2.294,00
 - 327 - 418 - 615 - 621 - 631

2294X2.05X1.07(GST)X1.15(freight)=

\$5600

YES	Code	Description	Additional price/ EURO
X			
	24	Arm supports fold-away, i.c.w. code 101/102	178,00
	29	Backrest, angle adjustable from seat onwards	180,00
	31	Backrest, angle reduction +/- 15° in 5° steps	180,00
	32	Seat width: 32 cm	0,00
	34	Seat width: 34 cm	0,00
	36	Seat width: 36 cm	0,00
	38	Seat width: 38 cm	0,00
	40	Seat width: 40 cm	0,00
	43	Seat width: 43 cm	0,00
	46	Seat width: 46 cm	0,00
	48	Seat width: 48 cm	0,00
	54	Footplate continuous 15 cm depth, 60 mm depth adjustable (in 4 steps), height-	218,00
	54/123	Footplate continuous 15 cm depth, 60 mm depth adjustable (in 4 steps), height-	236,00
	54/663	Footplate continuous 15 cm depth, 60 mm depth adjustable (in 4 steps), height-	290,00
	61	Driving wheels, wheel size 24 x 1"	0,00
	62	Wheel size 25 x 1" (pair)	0,00
	81	Arm support Aluminium, flip-up, incl. Height adjustable arm-pads, short	314,00
	84	Arm support Aluminium, flip-up, incl. height adjustable arm-pads, long	328,00
	87	Leg support 10 cm depth, 15 mm depth-, height- and angle adjustable, LLL 38-	236,00
	87/97	Continuous footplate, foldable, 10 cm width	366,00
	87/123	Leg support 10 cm depth, 15 mm depth-, height- and angle adjustable, LLL 30-	254,00
	87/663	Leg support 10 cm depth, 15 mm depth-, height- and angle adjustable, LLL 42-	308,00
	98	Wheel cover, 5 cm wide, i.c.w. Alu-Armrests Code 100, 102, 748, 749, 593	54,00
	100	Arm support Aluminium, reinforced, in frame colour, firmly screwed-on, with	154,00
	101	Arm support firmly screwed-on, with integrated clothes guard, black	148,00
	102	Arm support Aluminium, reinforced, in frame colour, firmly attached, with	224,00
	130	Colour: Titanium grey	0,00
	132	Colour: Red-metallic	0,00
	133	Colour: Green-metallic	0,00
	134	Frame Colour: Ocean blue	0,00
	140	Pushing handle short, telescopic, 10 cm height adjustable	170,00

YES	Code	Description	Additional
X			price/ EURO
	141	Push handle, foldable (BH 30 cm, not i.c.w. Code 502)	136,00
	143/900	Steering wheels 150 mm, puncture safe	0,00
	151	Steering wheel 100 solid rubber	0,00
	166	Handrim cover silicon, black, smooth	88,00
	168	Tetraplegic-silicon handrim 24" -Optigrip-	256,00
	173	Steering fork, short, for Code 327, incl. steering head, adjustable 46-52 cm	0,00
	183	Frame colour : Brilliant silver-metallic	0,00
	184	Frame colour : Bluestar	0,00
	187	Frame colour : Greystar	0,00
	188	Frame colour : Blackstar	0,00
	189	Frame colour : Violetstar	0,00
	203	Frame colour: Dark violet	0,00
	206	Frame colour: Matt black	0,00
	230	Frame colour: burgundy-red	0,00
	232	Cover colour: blue	0,00
	235	Cover colour: Deep black (for seat cushion)	0,00
	287	Cover colour: deep black (Nylon)	0,00
	292	Backrest, breathable open-neshed fabric, black	79,00
	300	Without steering wheel height adapter, FSH 52 cm	0,00
	301	Steering wheel height adapter, FSH max.55 cm	0,00
	305	Pressure Brake "Xchange" for driver	0,00
	309	Foot pedal tube	60,00
	312	Leg support continous, 10 cm, angle- und depth adjustable, (LLL 38-48 cm)	auf Anfrage
	312/663	Leg support continous, 10 cm, angle- und depth adjustable, (LLL 42-48 cm)	auf Anfrage
	313	Leg support continous, 15 cm, angle- und height adjustable (LLL 38-48 cm)	236,00
	313/663	Leg support continuous, 15 cm, height- and angle adjustable, flip-up and swivel away, LLL 38-48 cm	308,00
	325	Steering wheels 125 mm, solid rubber tyred	0,00
	327	Steering wheels 142mm soft solid rubber (i.c.w. Code 420 or 173)	0,00
	347	Compact frame X1 Speedster	0,00
	350	Seat width: 50 cm	328,00
	352	Frame length Standard 560 mm, Seat depth 36-42 cm	0,00
	353	Frame length long 600 mm, Seat depth 44 - 46 cm	0,00
	359	Calf belt (quick fastening)	0,00
	366	Seat depth: 36 cm	0,00
	368	Seat depth: 38 cm	0,00
	418	Steering wheel fork short, incl. steering head, adjustable 47-52 cm	0,00
	419	Steering wheel fork long incl. castor stem SH 49-55 cm	0,00
	420	Steering wheel fork long incl. steering head, adjustable 50-54 cm	0,00
	423	Wheel camber 0,5° factory-made	0,00
	445	Handrim titan with 6 holder (for active wheel)	256,00
	479	Hand and spoke guard, Decor "Xchange", only 24"	178,00
	486	Pneumatic tyres (standard profile, rear wheels)	0,00
	488	Ultra-light running tyres (driving wheels)	38,00
	489	Hand and spoke guard transparent	178,00
	493	Handrim stainless steel (pair)	0,00
	494	Handrim aluminium (pair)	0,00
	496	High pressure - tyres 8 bar, 24/25 x 1"	0,00
	502	Pushing handle, 20 cm height adjustable	148,00
	560	Seat cushion, 3 cm thick (textile)	80,00
	561	Seat cushion, 6 cm thick (textile)	80,00
	583	Arm support ALUMINIUM, reinforced, in frame colour, firmly screwed-on, with	314,00
	595	Butterfly brake for user	0,00
	615	Seat depth: 40 cm	0,00

YES	Code	Description	Additional
X			price/ EURO
	618	Seat depth: 44 cm	0,00
	619	Seat depth: 46 cm	0,00
	621	Backrest without lumbar support	0,00
	622	Backrest with lumbar support	0,00
	624	Spare code, select for seat height 39-43 cm	0,00
	626	Backrest height adjustable, 30-32,5 cm	0,00
	631	Backrest height adjustable: 35-37,5 cm	0,00
	632	Backrest height adjustable: 40-42,5 cm	0,00
	633	Backrest height adjustable: 45-47,5 cm	0,00
	651	Detachable driving wheels with quick release axle (pair) 24"	0,00
	664	Footboard locking	18,00
	679	Fastening strap	12,00
	681	Brake lever extension, swing-away	58,00
	691	Anti-tip roll, insert, swivel away (not i.c.w. Code 746)	192,00
	692	Handrim distance 15 mm	0,00
	693	Handrim distance 25 mm	0,00
	695	Active wheel incl. quick-release axle, radial spokes, high-flange hub, hollow	68,00
	696	Black wall tyres, Standard	0,00
	697	Tetraplegia device for silverdesign wheel (available from april 2007)	220,00
	712	Anti-tip roll, left mounted	96,00
	713	Anti-tip roll (right mounted)	96,00
	746	shunting castors	246,00
	748	Arm support ALUMINIUM, reinforced, in frame colour, firmly screwed-on, with	314,00
	749	Arm support ALUMINIUM, reinforced, in frame colour, firmly screwed-on, with	314,00
	784	Seat depth 42cm	0,00
	820	Speedster option backrest down-foldable (available from april 2007)	250,00
	823	Heel cup (pair)	0,00
	833	Safety belt	36,00
	878	Wheelbase extension 6,5 cm welded at rear frame	135,00
	883/	IRC Indoor tyres BLUE for 24" or 25"	
	854		145,00
	883/	IRC Indoor tyres RED for 24" or 25"	
	855		145,00
	895	Wheel width 1"	0,00
	896	Wheel width 1 3/8"	0,00
	897	Marathon Plus Evolution tyres, puncture safe 24"/25"	19,00
	899	Speedster option: wheel base extension Xtreme +15cm measured at backrest position, with second axle, detachable, restricted possibility to brake, incl. code 878. Available as of April 2007	280,00
	900	Puncture proof tyres for steering wheels	0,00
	902	Puncture proof PU - tyres drive wheels	0,00
	904	Puncture proof tyres for drive wheels (FLEXEL) 24"	166,00
	907	Negative camber 3° -factory made-	0,00
	911	Special air pump for wheelchairs with indoor sports tyres	92,00
	951	Accessory net	7,50
	962	Universal- rucksack with outer pocket	54,00
	970	Walking stick holder	36,00
	988	Leg protection cover (waterproof)	110,00
	4915	Frame colour: gold-special	120,00
	4916	Fram colour: titanium shine	0,00
	4930	Foot pedal tube, left assembly	30,00
	4931	Foot pedal tube, right assembly	30,00
	1054966	Occupant restraint system	285,00
	711	shop assembly	on request

YES	Code	Description	Additional
X			price/ EURO



x-trends, the trend items fitting for your



	1053166	Six-Star driving wheel for increased running, silence and stability. Incl. Aero rim, 8 bar high-pressure tyre and stainless steel handrim	188,00
	1052183	Six-Star driving wheel for increased running, silence and stability. Incl. Aero rim, 8 bar high-pressure tyre and stainless steel handrim	188,00
	1055579	Six-Star driving wheel for increased running, silence and stability. Incl. Aero rim, 8 bar high-pressure tyre and Tetra handrim	223,00
	1052184	CROSS driving wheel, incl. Alu rim (12 mm quick release axle not included) and cross country tyres 1.9/2.0, 24", handrim stainless steel,	130,00
	1061419	"Spinnergy-wheel", light, sturdy, reliable and attractive, 24" with Tetra-handrim	398,00
	1061420	"Spinnergy-wheel", light, sturdy, reliable and attractive, 24" with titan-handrim	398,00
	1061428	"Spinnergy-wheel", 25" with Tetra-handrim, incl. Marathon Plus Evolution tyres, 1 piece	398,00
	1061429	"Spinnergy-wheel", 25" with titan-handrim, incl. Marathon Plus Evolution tyres, 1 piece	398,00
	1056842	"Ultra-Light", 24", with Alu-handrim, tyres, black stainless steel spokes (12 mm quick-release axle not included)	188,00
	1061427	"Black-Light" Ultra light rim, 24", with Alu handrim, no flat, abrasion resistant and light (quick-release axle not incl.)	188,00
	261519300	X-Treme Rucksack, in cool aluminium optic	29,00
	1052185	Lightwheels 80 mm (1 piece)	24,00
	1052186	Lightwheels 100 mm (1 piece)	24,00
	1052188	Steering wheel with Aluminium-rim, solid rubber 100 mm	56,00
	1052189	Steering wheel with Aluminium-rim, solid rubber 125 mm	56,00
	261519600	Black Wheel-Case	70,00
	261599900	X-Chrono brand watch	100,00
	261530500	Silver-Pen	16,00
	261530900	"X-TREND" X-Beam Gyrotwister	34,50
	261531003	"X-TREND" LED-Light "X-CHANGE", white	8,00
	261531004	"X-TREND" LED-Light "X-CHANGE", red	8,00