

# NETTI 1

A light all-round fixed frame wheelchair.



Netti I is an all-round chair with fixed frame. This means that you achieve a stable and comfortable seat unit with excellent manoeuvring qualities.

Netti I can be adjusted to a seat height of 31 cm without cushion. This gives adequate space between the ground and the seat plate, offering freedom of movement for the lower legs should the users wish to propel the chair using their feet.

Backrest with recline function comes as an accessory. This is a good alternative for users who manoeuvre the chair with one leg and need a negative angle on the backrest to "lock" the pelvis and isolate the muscle power being used to pull the chair forward. The release handle can be mounted on either of the armrests.

Max users weight 125 kg.

For seat widths 55 and 60 cm the Netti I is tested and approved for users up to 150 kg.

## ALU REHAB AS

Fabrikkeveien 22 Kvål  
NO-4323 Sandnes, NORWAY  
Phone: +47 51 68 18 88  
Fax: +47 51 68 18 89

## NETTI I COMES WITH THE FOLLOWING EQUIPMENT AS STANDARD

<ul style="list-style-type: none"> <li>→ <b>BACK</b></li> <li>→ <b>ARMRESTS</b></li> <li>→ <b>FOOT SUPPORT</b></li> <li>→ <b>BRAKES</b></li> <li>→ <b>MAIN WHEELS</b></li> <li>→ <b>PUSH RIMS</b></li> <li>→ <b>FRONT CASTORS</b></li> <li>→ <b>ANTI TIP</b></li> </ul>	<p>Standard back with Velcro. The angle of the backrest can be adjusted in 3 fixed positions of 107°, 97° and 83°. The backrest can be folded down for transport.</p> <p>With skirt guard, height adjustable and removable. The armrest pads can be adjusted.</p> <p>Fixed angle 75°.</p> <p>The brake is mounted on the top frame tube (adjustable). 22" with quick release (In seat width 55 cm and 60 cm it comes with 24")</p> <p>Black plastic.</p> <p>6" with quick release (In seat width 55 cm and 60 cm it comes with 7")</p> <p>Height adjustable and length adjustable. Can be turned around and used as step climber.</p>
-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------



The back is foldable and makes transport of the chair very simple.

The configuration of the wheelchairs differs from market to market. For other accessory: See pages 34-42.

## TECHNICAL SPECIFICATIONS:

CHASSIS SIZE	FUNCTIONAL SEAT WIDTH	SEAT DEPTH* Standard (Adjustable in depth)	SEAT HEIGHT**		BACK HEIGHT (From seat plate to top Velcro)	WEIGHT KG
			Standard (Without change of parts)	(With change of parts; 20"-24")		
35	37,5	43(39-46)	39(40-50)	31-47	48	21,3
38	40,5	43(39-46)	39(40-50)	31-47	48	21,4
40	42,5	43(39-46)	39(40-50)	31-47	48	21,5
43	45,5	43(39-46)	39(40-50)	31-47	48	22,0
45	47,5	48(45-51)	39(40-50)	31-47	48	22,1
50	52,5	48(45-51)	39(40-50)	31-47	48	22,3
55	57,5	48(45-51)	41(40-50)	31-47	48	25,2
60	62,5	48(45-51)	41(40-50)	31-47	48	25,5

\* The seat depth is measured without back cushion from front of the seat plate to the back tube placed in standard position.

\*\* The seat height is measured from floor to top front seat plate.



Netti I offers adequate space for the feet.



The seat height on Netti I can be adjusted to a low position.

