

ALU REHAB ATOS IV

- Light weight, heavy-duty aluminium frame with 140Kg load capacity.
- Flip-up armrests. Detachable & height adjustable legrests.
- Flat-free 24" rear tyre, heavy-duty 8" front castors.
- 44cm seat-depth, 49.5cm seat-height.



Alu Rehab AS (Norway), has developed wheelchairs for the international market since the start in 1989. The wheelchairs are being developed in cooperation with users, therapists, nursing personnel, designers and producers. The wheelchairs are being developed from the thought that the user should get as much comfort as possible, combined with good manoeuvring qualities. A common denominator for the chairs is: Norwegian design based on ergonomics, safety and good functionality. Most of the accessories fit all the models

<http://www.alurehab.no>

Frame Width	40cm (15.7")	44cm (17.3")	48cm (18.9")	51cm (20.1")	54cm(21.2")
Weight (Kg)	16.9	17.6	18.3	19.7	19.5

Rehab Mart