



Atos III

The Atos wheelchair is designed to be a tool for improvement of life quality for disabled persons. It shall be a product both user and surroundings can trust.

The Atos III is an allround, good quality wheelchair made for both in-door and out-door use.

It is designed both for self-propelling and for use with assistance using the push handles.

It is equipped with height adjustable, swingable and removable foot supports.

Possibility for swinging the armsupports backward, gives good opportunities for side transfer and for getting the wheelchair under tables enabling activity.

The back upholstery is adjustable for individual fitting. This enhance good seat position.

The wheelchair is easily folded.

Max user weight is 140 kg



Atos III folded

ATOS III COMES WITH THE FOLLOWING EQUIPMENT AS STANDARD:

BACK	Back with velcro straps. The back has upholstery over the back tubes. It has a pocket on the back.
ARM SUPPORTS	The arm support has fixed height and is swingable.
FOOT SUPPORTS	Fixed angle 75°. Swingable and removable. Footplates can be rotated and they are adjustable in height.
BRAKES	The brake is mounted on the top frame tube.
MAIN WHEELS	24" main wheels fixed with screw. Puncture proof.
PUSH RIMS	Plastic
FRONT CASTORS	7" front castors. Puncture proof.

OPTIONS;

Seating system and foot supports from the Netti program

TECHNICAL SPECIFICATIONS:

Chassis size	Seat depth*	Seat height**	Back height***	Total width	Transport width	Weight Kg
35	44	49,5	41	52,5	31	16,9
40	44	49,5	41	57,5	31	17,6
45	44	49,5	41	62,5	31	18,3
50	44	49,5	41	67,5	31	19,1
55	44	49,5	41	72,5	31	19,9

* Measured from front of seat to backrest.

** Measured from the floor to top of seat in the rear.

*** Measured from seat cover to top velcro.