

KARMA KM-8020 / KM-8520

Light Weight Detachable Wheelchair

Suitable for: heavy users, spinal injury & drivers



- Ultra-lightweight aircraft-grade aluminum alloy frame
- Standard seat width: 14", 16" or 18"
- Open-meshed, water permeable upholstery.
- Weight capacity: 130kg
- flat-free rear wheels
- Flip-back or detachable armrests and swing-away footrests
- Foldable backrest
- Seat height adjustable
- Weight: 14.6kg (KM-8520 16")

KM-8520F24



KM-8020F24



KM-8520F12

Options:











- ⊙ Attendant cable brake
- ⊙ 14" rear wheels
- ⊙ Detachable and washable cushion
- ⊙ Quick-release rear wheels
- ⊙ Height Adjustable Armrest
- ⊙ Plastic spokes
- ⊙ Bariatric (Extra-Large) :
KM-8020X;
Seat width 20" / 22" / 24".
Weight capacity 160kg
- ⊙ PVC (water-proof) upholstery

Specifications

	KM-8020	KM-8520
Seat Width	46cm(18") / 51cm(20") / 56cm(22") / 61cm(24")	35cm(14") / 41cm(16") / 46cm(18")
Front Castor / Rear Wheels	8" / 24"	6" / 22"
Armrests	Detachable	Flip-up
Rear wheels	Mag wheels	Spoke wheels
Armrests	Detachable	Flip-up
Brakes	Parking brakes	Parking and attendant cable brakes
Backrest	Rigid	Foldable

*Max. User Weight: 130kg

Unit: cm

				
40.5/46	43/43	62.5/66	35/35	23/23
				
104/104	50/50	39/39	86/86	15.1kg/15.8kg

* 160Kg for KM-8020-20" / 22" / 24"

± 1cm ± 0.5kg

Rehab Mart
One Stop Home Care
www.RehabMart.sg

Rehab Mart (Pte) Ltd 23 Tagore Lane #04-19 Tagore 23 Warehouse Singapore 787601
Tel: 62500555 Fax:62507576 Email: mail@rehabmart.sg
Retail Outlets (Rehab Mart) **Balestier** : 12 Boon Teck Rd, #01-02 Shaw Corner T: 63524400 **Chinatown** 100 Eu Tong Sen St #02-34 Pearls Ctr T: 63247077 **Next to Paya Lebar MRT**: 10 Eunos Rd 8 #B1-02 Sing Post Ctr T:64551551 **Upp Thomson**: 23 Tagore Lane #04-19 Tagore 23 Warehouse T:62500555