

Karma KM-AA20

Wheelchair for Airport & Airplane Aisle



- Ultra-lightweight aircraft-grade aluminum alloy frame
- Weight: 11.5kg (w/o rear wheels); 16.2kg (with rear wheels)
- Weight capacity: 100kg
- Swing-away armrests for convenient wheelchair entrance and exit
- Sturdy non-foldable frame
- Foot-operated quick-release rear wheels for easy airplane boarding and unboarding
- Wooden armrests and security chest and waist belts
- Can be customized according to airplane aisle width

KM-AA20 Specifications

Unit: cm

*Size: 14" *F/R Wheels: 6"/24" (7") *Max. User Weight: 100kg

39	42	58 (54)	-	19
99 (87)	48.5	41.5	90.5 (88.5)	16.2kg (11.5kg)
				±1cm ±0.5kg



Rehab Mart