

Karma KM-1000 / KM-1500

Standard Light Weight Aluminium Wheelchair

For general users, healthcare centers, and hospitals



■ Ultra-lightweight aircraft-grade aluminum alloy champagne-coloured frame.

■ Upholstery :

- KM-1000 - Water-proof PVC
- KM-1500 - Open-mesh, breathable.

■ Seat width:

- KM-1000 - 18"
- KM-1500 - 16" & 18"

■ Foldable frame

■ Weight capacity: 100kg

■ 8" casters and 24" flat-free rear wheels

■ Fixed armrests and footrests

■ Backrest :

- KM-1000 - Rigid
- KM-1500 - Foldable

■ Weight: 13.5kg

Options with extra costs:

- ◎ Attendant cable brake
- ◎ Detachable and washable cushion



KM-1500 F14 (push-chair)

Specifications :

*Size: 16"/18"

*Wheels: 8"/24"

*Max. User Weight: 100kg

39.5cm/45.5cm	41cm/41cm	62cm/68cm	31.5cm/31.5cm	19cm/19cm
101cm/101cm	50cm/50cm	37cm/37cm	87cm/87cm	13.5kg/13.6kg

± 1cm

± 0.5kg

Rehab Mart