

CUSTOMIZE!



Amigo offers choices. YOU CHOOSE, the color, wheels, seat and accessories for the Amigo based on your needs.

Click [here](#) to read what one happy owner wrote about his customized Amigo.

Click on any selection to view an enlarged photo with description, part number and MSRP. Give our Sales Department a call for any further assistance at 800-248-9130. They'll be pleased to assist you!

Customized Options		
	Adapt-A-Seat Frame Custom Dimension (Part Number: 396000) MSRP: \$600.00	<ul style="list-style-type: none"> Accommodates High-End Rehab Seating Systems Up To 24" Wide Steel Framework For Custom Seating
	Adapt-A-Seat Pan Custom Dimension (Part Number: 396001) MSRP: \$140.00	<ul style="list-style-type: none"> Accommodates High-End Rehab Seating Systems Up To 24" Wide Steel Framework For Custom Seating
	Arm Assembly, Height Adjustable (pair) (Part Number: 394060) MSRP: \$300.00	<ul style="list-style-type: none"> Rest Your Arms At Just The Right Height Adjusts From 1" Lower Than Standard to 3" Higher Than Standard
	Basket, Large (Part Number: 398073) MSRP: \$200.00	<ul style="list-style-type: none"> Rear Mount Removable Convenient Size & Location 16"w x 9"d x 14"h
	Basket, Small (Part Number: 398072) MSRP: \$100.00	<ul style="list-style-type: none"> Rear Mount/Side Mount Front Mount Alternative 13.5"w x 10"d x 10"h
	Calf Pad, Single (Part Number: 398005) MSRP: \$350.00	<ul style="list-style-type: none"> Calf Pad And Support Assembly (Single or Pair)

	<p>Flag, Safety and Holder (Part Number: 398074)</p> <p>MSRP: \$65.00</p>	<ul style="list-style-type: none"> • Rider More Visible & More Easily Located at Longer Distances • 70" High Pole
	<p>Foot Plate, Adjustable (Part Number: 399999)</p> <p>Call for Quote</p>	<ul style="list-style-type: none"> • Comfort Adjustment-Vertically & Horizontally • Seat Bottom Attachment - Barber Style • Foot Plate Moves With You
	<p>Foot Retainer w/o Strap (Part Number: 398023)</p> <p>MSRP: \$55.00</p>	<ul style="list-style-type: none"> • Keeps Feet Resting On Platform
	<p>Foot Retainer w/Strap (Part Number: 398022)</p> <p>MSRP: \$85.00</p>	<ul style="list-style-type: none"> • Keeps Feet Still • Positioned Safely & Correctly
	<p>Horn, Automatic Reverse Alert (Part Number: 398067)</p> <p>MSRP: \$90.00</p>	<ul style="list-style-type: none"> • Sounds Reverse Movement
	<p>Leg Retainer for Flip Up Arms (pair) (Part Number: 398002)</p> <p>MSRP: \$240.00</p>	<ul style="list-style-type: none"> • Upper Thighs and Hips Stay In Tight/Secure
	<p>Oxygen Holder, Dual Tank (Part Number: 399999)</p> <p>Call for Quote</p>	<ul style="list-style-type: none"> • Two E-Cylinders of O2 • Longer Distances Coverage
	<p>Pin, Fixed Arm (Part Number: 399999)</p> <p>Call for Quote</p>	<ul style="list-style-type: none"> • Eliminates Flip-Up Capability
	<p>Power Seat Switch, Remote (Part Number: 395010)</p> <p>MSRP: \$320.00</p>	<ul style="list-style-type: none"> • Used on Rear Drive and RT Express Amigos • Toggle Switch Used w/HI-LO Speed Knob • Assists w/Less Manual Dexterity
	<p>Seat Cover, (Back & Bottom Cushions) Hand-Crafted Leather (Part Number: 396015)</p> <p>MSRP: \$400.00</p>	<ul style="list-style-type: none"> • Upholstery Grade • Made From Genuine Leather Hides • Backed w/Quilting Foam For Finished Look • Professional Grade High Density Foam Padding
	<p>Seat Post Assembly, Lower (Smooth Ride) (Part Number: 397006)</p> <p>MSRP: \$350.00</p>	<ul style="list-style-type: none"> • Absorbs Jarring Rougher Terrain • Capacity 175lbs.+
	<p>Seat Post, (Shorter) Power Seat Lift (Part Number: 397015)</p> <p>MSRP: \$1260.00</p>	<ul style="list-style-type: none"> • For 24V RD System • Available At 1" Lower Than Standard

	<p>Seat Post, (Shorter) Power Seat Lift (Part Number: 397020)</p> <p>MSRP: \$1260.00</p>	<ul style="list-style-type: none"> • For 12V FD System • Available At 1" Lower Than Standard
	<p>Seat Post, Power Seat Lift (Part Number: 397011)</p> <p>MSRP: \$1700.00</p>	<ul style="list-style-type: none"> • 2-Stage 24V 7.75" • Maximum Vertical Travel • Seat Height 22" to 29-3/4" • 175lb. Weight Capacity
	<p>Seat Post, Power Seat Lift (Part Number: 397012)</p> <p>MSRP: \$1700.00</p>	<ul style="list-style-type: none"> • 2-Stage 12V 5.75" • Starting At 1" Lower Than Standard Seat Post • Seat Height 21" to 26-3/4" • 250lb. Weight Capacity • Greater Vertical Travel
	<p>Seat Post, Power Seat Lift (Part Number: 397013)</p> <p>MSRP: \$1700.00</p>	<ul style="list-style-type: none"> • 2-Stage 24V 5.75" • Starting At 1" Lower Than Standard Seat Post • Seat Height 21" to 26-3/4" • 250lb. Weight Capacity • Greater Vertical Travel
	<p>Seat, Custom Dimensions (Part Number: 396005)</p> <p>MSRP: \$800.00</p>	<ul style="list-style-type: none"> • Up to 24" Wide • Cloth & Vinyl • Wide Seat Comfort • Metal Seat w/Spring Bottom Cushion
	<p>Step, Platform (mini) (Part Number: 393000)</p> <p>MSRP: \$350.00</p>	<ul style="list-style-type: none"> • For Shorter Legs & Feet • Height Selected Step As Required
	<p>Switch, Rocker On/Off (Part Number: 398068)</p> <p>MSRP: \$30.00</p>	<ul style="list-style-type: none"> • Eliminates Need For Key
	<p>Throttle Extensions (pair) (Part Number: 392022)</p> <p>MSRP: \$85.00</p>	<ul style="list-style-type: none"> • Wider Tabs for Easier Palm Use • Complements Lesser Manual Dexterity
	<p>V-Trak Backrest System (Part Number: 10753, 10754, 10755 & 10757, 10756 & 10757)</p> <p>MSRP: \$2670.00</p>	<ul style="list-style-type: none"> • Back 17" w x 23" h, BM L 17 • Headrest HS AW • Right Side Lateral Ts2 RH w/Ts2 MP Pad • Left Side Lateral Ts2 LH w/Ts2 MP Pad
	<p>Wheels, Stability/Steel Platform (pair) (Part Number: 398075)</p> <p>MSRP: \$100.00</p>	<ul style="list-style-type: none"> • Additional Anti-Tip Prevention

Specialized Customization Selections

Samplings Requested By Customers:

- Barber-Style Adjustable Foot Rest
- Seat Extension Pads
- Scoliosis Pads
- Leg Rests
- Trailer With Basket
- 4" Forward Arm Pads
- 1-3/4" Shorter Handle
- 12V Accessory Outlet Installed
- Rear View Mirror
- Quad Cane Holder
- Platform Dimensions (length only)
- Automotive Color Selections (specific PPG number required)