



# STRUGGLING TO FIND A HIGH QUALITY MOBILITY PRODUCT AT A LOW PRICE?

## The Amigo MC Is Your Best Value Solution

- Industry-Recognized Quality at Low Pricing
- Robust Weight Capacity
- Durable Rear Drive Performance
- Two Model Selections



**NEW!**  
**X**FACTOR

**MCX - #690000B**  
**ONLY \$1995 MSRP\***  
*(complete w/batteries)*



**MC - #310000B**  
**ONLY \$1995 MSRP\***  
*(complete w/batteries)*

# MC & MCX



**OTHER CUSTOMIZED MOBILITY SOLUTIONS**



**Upgraded MC Shown**  
Options: Premier I seat, 10.5" Trimline  
Tires & Cinnamon Red Handle



**Upgraded MCX Shown**  
Options: Premier I seat, 9" Trimline  
Tires & Scuba Blue Handle

## Features & Benefits

	<b>MC</b>	<b>MCX</b>
Solid Steel Platform and Frame-Long Lasting Durability	✓	✓
On-board Diagnostics-Heads-Up On Service Issues	✓	✓
Fully Rotational Seat of 360 Degrees-Position Exact Setting of Choice	✓	✓
Robust Weight Capacity-Generous Comfort	350 lbs	250 lbs
Infinitely Adjustable Seat Lock-Lock Seat In Select Position	✓	✓
Unique Oval Handle-Easy Grip/Perfect for Limited Hand Usage	✓	✓
Infinitely Adjustable Black Painted Handle-Comfortable Reach	26"(66.04cm)	25" (63.5cm)
Tight Turning Radius-Easy Maneuverability	34"(86.4cm)	33"(83.8cm) <b>XFACTOR</b>
Color Selection	Cinnamon Red	Scuba Blue
Battery/Travel Distance	(2) 17amp Batteries 18-23 Miles (29-35km)	(2) 17amp Batteries 22 Miles (35.4km)
Range of Speed	0 - 5.9mph (9.5kph)	0-5.0mph (8.0kph)
Dynamic and regenerative brakes	✓	✓
Flat-free indoor/outdoor tires	✓	✓
Standard Universal Front Basket	✓	✓
Amigo limited lifetime/3 year warranty	✓	✓
HI/LO Speed Controls	✓	✓
2 AMP-24V Off-board Universal Charger w/AC Cord	✓	✓
Flip-up, Adjustable, Removable Arms	✓	✓
Length	43.75" (111.1cm)	41.5" (105.4cm) <b>XFACTOR</b>
Rear Wheel Width	23.50" (59.7cm)	20.5" (52.1cm) <b>XFACTOR</b>

### Amigo Lift Options: Take It With You

Amigo Lift Options easily transport Amigo products in a majority of vehicles.  
Contact your local dealer for additional information and installation requirements.

*\*Prices subject to change without notice (US Dollars).*

Your Local Representative:



**AMIGO MOBILITY INTL, INC.**  
(800) 248-9130 • (989) 777-0910  
[www.myamigo.com](http://www.myamigo.com)

© 2006 Amigo Mobility International, Inc. Amigo® is a trademark of Amigo Mobility International, Inc.