

## KP-45.3TR

The model designed for active outdoor use.

The KP-45.3 TR is a stylish, high performance powered wheelchair offering excellent maneuverability. Equipped with a powerful, high-efficiency 400W motor, ergonomically designed seat and armrest enabling both tilt and recline adjustment for optimal comfort, and a wide variety of customizing options, this is the chair you've been waiting for.



KP-45.3 TR equips with multi-functional adjustable seat, the backrest is capable of reclining up to 122°, offering pressure-relieving comfort at will.

■ Fully adjustable seat & armrest system for great comfort on your demand



Forward / backward adjustment



Height adjustment



Fully detachable



Width adjustment

■ Colors



### Features

- 400W high performance motor
- Built-in light set design
- Suspension anti-tippers

- Multi-dimension adjustable seat & armrest
- Recline backrest
- Built-in charger socket protecting the controller



Multi-functional controller



Freewheel hub



Rear wheel fender



Suspension anti-tippers



Built-in lighting system, improve driving safety at night

### Optional accessories

- 75 Ah battery
- High-speed tires
- Parallel swing-away joystick bracket

- L-shaped hand bar
- Curb climber (please notify us upon placing order)
- Drink holder



Parallel swing-away joystick bracket



L-shape hand bar



Push bar

Power tilt and recline allows user to relieve pressure while adjusting to the most comfortable sitting position, effectively spreads weight to a greater area and reduce the risk of decubitus ulcers.



Traditional powered Wheelchair - High risk of causing decubitus ulcers because of pressure points concentrating at the bottom.

Powered tilt & recline wheelchair - Power controlled tilting and reclining enable seat angle adjustment to optimize comfort, spreading pressure of bottom and back at ease.

■ The attributes mentioned above are only suitable for KP-45.3 TR & KP-45.3 T

### KP-45.3T

KP-45.3 T powered backrest enables tilting and reclining simultaneously, greatly relieves pressure and prevents decubitus ulcers from occurring.



Power tilt up to 38°



Manual recline up to 122°

### KP-45.3

Quick seat detachment system, convenient for transport & easy maintenance.



Quick seat detachment design (KP-45.3 Only)



Manual recline up to 122°








FDA MDA  CE ISO 9001

104-40001 Rev.02.AUG.2008

Discover the path and make it yours.

**KARMA**<sup>TM</sup>  
mobility solutions

															
18"/20"	122cm	65cm	127cm	Front Wheel/Adjustable	400Wx2	50Ahx2	Shark68A	90cm	8"	13km	40-50km	Aluminum Alloy	135kg	107kg	78kg

The above are estimated data under 20°C (68°F), an 80kg (176lbs) user, brand new fully charged batteries and a constant driving speed of 6km/hr (3.7mph). Distance, climbing ability and operating time can vary according to condition.

± 1cm ±0.5kg



**KP-45.3 Series**

The model designed for active outdoor use

No. 2363 Sec. 2 Da-Shiue Rd.,  
Min-Hsiung Shiang, Chia-Yi 621, Taiwan  
Tel: 886-5-206-6688 Fax: 886-5-206-7788  
E-mail: globalsales@karma.com.tw

**KARMA**<sup>TM</sup>  
mobility solutions

KARMA MEDICAL PRODUCTS CO., LTD.

www.KarmaMedical.com