

Power Wheelchairs

This foldable wheelchair can be easily carried and stored in any vehicle for active and independent lifestyles. A single person can easily fold and store the power chair without the use of any lifting device. Ideal for indoor use with its compact size and excellent maneuverability

Features:

- Easily foldable battery bracket allows improved portability
- Lightweight and foldable for easy handling, storage and transportation
- Programmable controller with a wide range of selections to suit the needs and demands of different users
- Rugged heavy-duty casters
- Auxiliary brakes offer extra safety
- Adjustable height flip-up armrests, swing away footrests and optional legrests enable close access to objects
- 23" turning radius allows the chair to be used in very confined spaces
- Optional PG Pilot Controller
- Limited lifetime warranty on frame

Options(inch):

Model	Seat		Arms		Back	Footrests		Overall dimensions				W/batteries Loading W.	
	W	D	W	D	H	W	D	W	D	W	H	W	
HE98000001	18	18	20	20	9-14	19	15	19	40	37	24.5	14	67Kg 130Kg
HE98000002	18	18	20	20	9-14	19	15	19	44	37	27.5	14	67Kg 130Kg

Technical Specs:

Model	HE98000001	HE98000002
Max. Speed	6.5kph/4 mph	8/5mph
Range	32km/20mi.	32km/20mi.
Turn radius	58cm/23"	58cm/23"
Motor	DC 24V, 180W, 3800rpm	DC 24V, 200W, 2400rpm
Controller	PG, VSI, 50A	PG, VSI, 50A
Battery	12V, U1(36AH) x2	12V, U1(36AH) x2
Charger	5 Amp off-board	5 Amp off-board
Gradient	12°	12°
Caster	8" (200 x 50) foam filled tire	8" (200 x 50) foam filled tire
Drive wheel	12 1/2" foam filled tire	24 x 1-3/8" foam filled tire
Brake	Intelligent, regenerative, Electromagnetic brakes	Intelligent, regenerative, Electromagnetic brakes



HE98000001



HE98000002