



# Product Catalogue

20100527



# C300

With the C300 you acquire a front-wheel drive Permobil with all you could desire in the way of a safe wheelchair, and that at a price which matches your wallet. A slightly smaller and more flexible chassis, perfect for indoor use, but still a strong and robust wheelchair with good outdoor performance. We are talking Scandinavian design combined with cutting end technology and a proven staying power.

The C300 comes with four different seating systems. The simpler systems provide comfortable, slightly contoured seats for clients who feel comfortable and stable with a low profile seat. For those who have high demands on the seat the Corpus seat will be the natural choice. A classic with few competitors in terms of ergonomics and adaptability. Read more about all our seating systems, at the Seating systems section.

Max. speed:	7 km/h
Range:	30-35 km
Obstacle clearance ability:	70 mm
Colors available:	red, blue, green, black metallic
Width:	620 mm
Length:	1080 mm

Seat tilt    Seat lift    Leg rest    Back rest

Turning radius 625 mm    Seat height 510-710 mm



## C300 Corpus

Ergonomically approved  
 Modules allow for changes to be made  
 Can be fitted with all accessories and additional functions  
 Rnet Electronics  
 Max user weight 120kg

Turning radius 625 mm    Seat height 510-710 mm



## C300 PS

A smart seat that is easy to adjust  
 A sliding, telescopic frame and expanding seat pan  
 VR2/Rnet Electronics  
 Max user weight 120kg

Turning radius 625 mm    Seat height 460-660 mm  
 510-710 mm



## C300 TS

Built-in UniTrack  
 Seamless adjustment of seat depth and arm rest height  
 Functional design  
 VR2 Electronics  
 Max user weight 136kg

Seat tilt    Seat lift    Leg rest    Back rest



Turning radius  
625 mm



Seat height  
490-690 mm



## C300 CS

Basic but comfortable seat with manual adjustment of backrest angle  
VR2 Electronics  
Max user weight 120kg

## Accessories chassis



### Lights

1821059, VR2  
1824875, R-net



### Cover

Strong nylon  
310497-99-0



### Permlock C3

It allows you to drive or ride while sitting in your Permobil.

### Locking device

1823300

### Locking pins

Mounts underneath your wheelchair  
1824229

### Battery 60 Ah

101969-99-0

### Solid wheel

103499-99-0

### Mud and Snow tire

104871-99-0

### Charger JK Medico 10A

1824895

### Bag Charger

1824898