



## Netti 4U comfort|CE

- Compact, low-weight and narrow total width
- Simple in use
- Excellent functional design in the ergonomic solutions

Netti 4U comfort|CE is a designer wheelchair – created to meet the demands of comfort and relaxation. The combination of ergonomic shaped seat unit and larger number of adjustment possibilities allows comfortable and individual configuration.

The arm support and foot support are easily removed and give good opportunities for side transfers.

Since the total width is narrow and the chair has good manoeuvring performance – the chair is especially good for rooms with limited space.

Low weight and easy folding makes the wheelchair easy to handle.

The following accessories are available for both Netti 4U comfort|CE and CES. Several variants of head support and hip belts, abduction block in small, medium and large. Anti tip with pedal and 12", 16" and 24" main wheel can also be fitted.

## Netti 4U comfort|CES

Netti 4U comfort|CES is identical with the Netti 4U comfort|CE but 50mm shorter. Some parts will therefore differ.



Netti 4U comfort|CES tilted backwards.



# Rehab Mart

Measures and weights – look at page 31

## Netti 4U comfort|CED

- Compact, low-weight and narrow total width
- Simple in use
- Excellent functional design in the ergonomic solutions

Netti 4U comfort|CED is a multi-functional high-comfort wheelchair with excellent ergonomic solutions, narrow total width and low total weight. The simple functional controls are marked red and simple adjustment facilities make this chair extremely user-friendly.

The sitting unit is individually adapted and provides unique stability, support and comfort.

The wheelchair is small (narrow total width) and compact and is therefore very suitable for rooms with limited space of narrow door openings.

The wheelchair's low weight and folding frame make it easy to transport.

Removeable or revolvable arm support allow for simpler and better usage of hoist and sideways transfer.



Netti 4U comfort|CED is available both with angle adjustable and universal foot support.

A wide choice of accessories is available for this chair.



Revolvable arm support allow simpler and better sideways transfer.



The unit gives unique stability, support and comfort.



The wheelchair is small and compact and is suitable for rooms with limited space.



Low weight only 29 kg complete, together with folding frame, make it easy to transport.

Measures and weights – look at page 31