



Motorized Standing Wheelchair LY-ESB-220

Specification:

Motor: 450Wx24Vx2pcs

Battery: 50AHx12Vx2pcs

Charger: 24V x 5AMP.

Max loading: 120kg

Front wheel : 8" x 2" PU castor

Brake: Magnetic brake

Travel range: up to 30km

Max speed: 8.5km/h

Anti-tilted side support

Anterior/posterior stability mechanism

Electric controlled standing mechanism

Actuator: 12/24V DC permanent magnetic motor

Model	Size	PACKING SIZE (L x W x H) (inch)								Front castor (inch)	Rear wheels (inch)	N.W(kgs.) (incl. battery)	
LY-ESB22LR	XL	50" X 31" X 31"	1210	740	990	540	480	500	148-287	540	8	12	98.4
LY-ESB22AR	L	49.5" X 29.5" X 31"	1210	700	950	540	455	460	148-287	520	8	12	95
LY-ESB22DR	M	48" X 27.5" X 31"	1180	660	910	525	425	420	148-287	450	8	12	91.5
LY-ESB22SR	S	48" X 26" X 31"	1180	620	900	525	400	380	148-287	440	8	12	89.5