



# Karma SME

## Standing Wheelchair

Manual wheelchair with power control for standing up

- Four-points touch ground in standing position providing extra safety
- Ultra-light weight--80% of the frame is made of ultra light aircraft-grade aluminum
- Adjustable seat depth from 44cm, 47cm
- Seat width: 40cm, 44cm
- Ergonomic flip back armrest with height adjustable range from 18cm to 26cm
- Foldable footrest with angle adjustment characteristics
- Angle-Adjustable knee pad to fit individual preference

## SME Standing Wheelchair

The ability to stand up strongly promotes user independence, increases health benefits and improves life imbalances. The user will feel comfortable and confident in a near-natural standing position.

### Benefits

Activates circulation system and improve cardio-pulmonary functio

Helps prevent decubitus sores

Improves bowel regularity and function











Helps bladder and intestinal motility

Reduces contracture and ankylosis

Helps in preventing atrophy of the leg muscles

### Specification

Unit: cm

Model	Width	F/R wheels				
SME	40/44		40/44	44/47	59/63	
						Max. User Weight
	101.5	51	42	96/143	28.6KG	120KG

**Rehab Mart**  
[www.RehabMart.com.sg](http://www.RehabMart.com.sg)  
 One Stop Home Care

**A division of Greentyle Pte Ltd**

**Office:** 12 Boon Teck Rd, #01-02 Shaw Corner, S'pore 329586

Tel: 62500555 Fax: 65-62507576 Email: [rehab@pacific.net.sg](mailto:rehab@pacific.net.sg)

**Retail: Balestier:** 12 Boon Teck Rd, #01-02 Shaw Corner T: 62500555

**Outram Park MRT:** 100 Eu Tong Sen St #02-34 Pearls Ctr T: 63247077

**Paya Lebar MRT:** 8 #B1-02 Sing Post Ctr T:6455 1551

**SGH:** Blk 7 L1 T:6221 1110